

INTRODUCING

SUREFUSER™ +

VARIABLE AND CONTINUOUS ELASTOMERIC INFUSION SYSTEM



A solution for ambulatory infusion therapies

For use with various chemotherapy drugs, thalassemia treatment, pain medication and antibiotics

fluorouracil / adruciladrucil (5 -fu) -most used drug

Clear and easy read-out of the therapy progression

blue progression line moves along the scale

Safe: back-check valve with integrated infusion line

reduces risk for leakage and contamination

Strong lightweight housing

protects the contents from direct heat and pressure that can influence the flow rate

Dual filter

removes air from the infusion line and prevents particles larger than 0.2 µm from passing through

Products registered and listed on SAOC, ICON, and MEDPRAX:

Reference code	Description	Max volume	Flow rate	Nappi code (Singles & Box)	Box Qty
SFS-1001DP-EJ	Surefuser home infusion pump 100 ml/24 hrs	100 ml	4,2 ml/hr	1102842001 1102842002	1 10
SFS-1002DP-EJ	Surefuser home infusion pump 100 ml/48 hrs	100 ml	2,1 ml/hr	1102843001 1102843002	1 10
SFS-215P-EJ	Surefuser home infusion pump 150 ml/75 hrs	150 ml	2 ml/hr	1102844001 1102844002	1 5
SFS-525P-EJ	Surefuser home infusion pump 250 ml/50 hrs	250 ml	5 ml/hr	1102845001 1102845002	1 5
SFS-530P-EJ	Surefuser home infusion pump 300 ml/60 hrs	300 ml	5 ml/hr	1102846001 1102846002	1 5

Call your Nipro representative for more information and training:
Tharine De Beer | 066 327 5703 | tharine.debeer@nipro-group.com

